

Product Datasheet

**TOPO II Alpha (phospho-Ser1106) antibody (PE)
(orb504838)**

Description

Rabbit polyclonal antibody against

conjugated to PE

Species/Host

Rabbit

Reactivity

Human

Conjugation

PE

Tested

ICC

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human TOPO II Alpha around the phosphorylation site of Ser1106

Target

TOP2A

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

174kDa

Uniprot ID
P11388
Dilution Range

IF: 1:50-200

Expiration Date

12 months from date of receipt.