

Product Datasheet

IKK gamma (phospho-Ser31) antibody (PE) [Out of stock] (orb504949)

Description

Rabbit polyclonal antibody against

conjugated to PE

Species/Host

Rabbit

Reactivity

Human

Conjugation

PE

Tested

FC

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human IKK gamma around the phosphorylation site of Ser31 EE(p-S)PL

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

46kDa

Uniprot ID
Q9Y6K9
Dilution Range

IF: 1:50-200

Expiration Date

12 months from date of receipt.