

Product Datasheet

Histone H3 (phospho-Thr11) antibody (PE) (orb504970)

Description

Histone H3 (phospho-Thr11) antibody

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PE
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Histone H3 around the phosphorylation site of Thr11
Target	Phospho-Histone H3 (Thr11)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	Histone H3.1
Uniprot ID	P68431
Dilution Range	Flow-Cyt=0.2µg /Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.