

Product Datasheet

c-Fos (phospho-Ser32) antibody (PE) (orb504977)

Description	c-Fos (phospho-Ser32) antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	PE
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human c-Fos around the phosphorylation site of Ser32
Target	Phospho-c-Fos (Ser32)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01100
Expiration Date	12 months from date of receipt.