

Product Datasheet

Cyclin E(phospho-Thr62) antibody (PE) (orb504990)

Description	Cyclin E(phospho-Thr62) antibody (PE)
Species/Host	Rabbit
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Cyclin E around the phosphorylation site of Thr77
Target	phospho-Cyclin E1 (Thr77)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P24864
Expiration Date	12 months from date of receipt.