

Product Datasheet

Androgen Receptor (phospho-Ser595) antibody (PE) (orb505044)

Description

Androgen Receptor (phospho-

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Androgen Receptor around the phosphorylation site of Ser597
Target	phospho-Androgen Receptor (Ser597)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	b3
Uniprot ID	P10275
Dilution Range	IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.