

Product Datasheet

**AKT (Phospho-Tyr473) antibody (PE) [Out of stock]
(orb505129)**

Description

Rabbit polyclonal antibody against

PE

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

PE

Tested

FC

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human AKT3 around the phosphorylation site of Tyr473 FS(p-Y)SA

Preservatives

0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

57kDa

Uniprot ID
P31749
Dilution Range

ICC: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.