

Product Datasheet

STAT3 (phospho-Tyr705) antibody (PE) [Out of stock] (orb505290)

Description

Mouse monoclonal antibody against

conjugated to PE

Species/Host

Mouse

Conjugation

PE

Tested

ICC, IF

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human STAT3 around the phosphorylation site of Tyr705

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

MW

88kDa

Uniprot ID

P42227, **P40763**, **P52631**

Dilution Range

ICC: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.