

## Product Datasheet

**BACE1 (phospho-Ser498) antibody (PE) [Out of stock] (orb505391)**

**Description**

Rabbit polyclonal antibody against

conjugated to PE

**Species/Host**

Rabbit

**Conjugation**

PE

**Tested**

ICC, IF

**Applications**
**Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human BACE1 around the phosphorylation site of Ser498

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Lyophilized or Liquid

**Concentration**

- 100 µg (in 100 µl): 1 mg/ml

**Storage**

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

51kDa

**Uniprot ID**
**P56817, P56819, P56818**
**Dilution Range**

ICC: 1:50-200, IF: 1:50-200

**Expiration Date**

12 months from date of receipt.