

Product Datasheet

Caspase 9 (phospho-Thr125) antibody (PE) (orb505607)

Description	Caspase 9 (phospho-Thr125) antibody (PE)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PE
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Caspase-9 around the phosphorylation site of Thr125
Target	Phospho-Caspase-9 (Thr125)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P55211
Dilution Range	Flow-Cyt=0.2µg /test,IF=1:100-500
Expiration Date	12 months from date of receipt.