

Product Datasheet

PDE2A antibody (Biotin) (orb50561)

Description	PDE2A antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human cGMP-dependent 3',5'-cyclic phosphodiesterase protein (15-226AA)
Target	PDE2A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O00408
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.