

## Product Datasheet

### p90RSK (phospho-Thr573) antibody (PE) [Out of stock] (orb505714)

|                            |                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | p90RSK (phospho-Thr573) antibody (PE) [Out of stock]                                                           |
| <b>Species/Host</b>        | Rabbit                                                                                                         |
| <b>Conjugation</b>         | PE                                                                                                             |
| <b>Tested Applications</b> | ELISA, IF, IHC-Fr, IHC-P, WB                                                                                   |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from human RPS6KA1 around the phosphorylation site of Thr573 |
| <b>Target</b>              | Phospho-RPS6KA1 (Thr573)                                                                                       |
| <b>Preservatives</b>       | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.                                          |
| <b>Form/Appearance</b>     | Liquid                                                                                                         |
| <b>Concentration</b>       | 1mg/ml                                                                                                         |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                               |
| <b>Note</b>                | For research use only                                                                                          |
| <b>Isotype</b>             | IgG                                                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                                                     |
| <b>Uniprot ID</b>          | <b>Q15418</b>                                                                                                  |
| <b>Dilution Range</b>      | WB: 1:500-2000, IHC-P: 1:100-500, IHC-F: 1:100-500, IF: 1:100-500, ELISA: 1:5000-10000                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                |