

Product Datasheet

PDK1 (phospho-Tyr9) antibody (PE) [Out of stock] (orb505921)

Description	Rabbit polyclonal antibody against PDK1 (phospho-Tyr9) conjugated to PE
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	PE
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PDK1 around the phosphorylation site of Tyr9 QL(p-Y)DA
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61kDa
Uniprot ID	O15530
Dilution Range	ICC: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.