


Product Datasheet

Tubulin beta 3 (phospho-Ser172) antibody (orb506110)

Description	Tubulin beta 3 (phospho-Ser172) antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human TUBB3 around the phosphorylation site of Ser172 (VP(p-S)PK)
Target	phospho-TUBB3 (Ser172)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50-55 kDa
Uniprot ID	Q13509
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Rat brain
using
phosph...