

Product Datasheet

Bad (phospho-Ser118) antibody (orb506123)

Description	Bad (phospho-Ser118) antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Bad around the phosphorylation site of Ser118 (RM(p-S)DE)
Target	Phospho-Bad (Ser118)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	22 kDa
Uniprot ID	Q92934
Expiration Date	12 months from date of receipt.