



Product Datasheet

beta Actin (phospho-Tyr53) antibody (orb506128)


Description	beta Actin (phospho-Tyr53) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human beta Actin around the phosphorylation site of Tyr53 (DS(p-Y)VG)
Target	phospho-beta Actin (Tyr53)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42 kDa
Uniprot ID	P60709
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1ug/Test, IHC-F=1:100-500, ICC=1:100, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow cytometry of
NIH/3T3.Primary
Antibo...



Western blot:Lane
1:HeLa (Human)
Cell;La...



Western
blot:Pancreas
(Mouse);Primary:An...