

Product Datasheet

Gli1 (phospho-S112 antibody (orb506130))

Description	Gli1 (phospho-S112 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated sSynthesised phosphopeptidepeptide derived from human Gli1 around the phosphorylation site of S112/ T115/ S116 (IN(p-S)RC(p-T)(p-S)PG)
Target	phospho-Gli1 (S112/ T115/ S116)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/1ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	118 kDa
Uniprot ID	P08151
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.