



Product Datasheet

PERK (phospho-Thr980) antibody (orb506147)


Description	PERK (phospho-Thr980) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from mouse PERK around the phosphorylation site of Thr980 (H(p-T)GQ)
Target	Phospho-PERK (Thr980)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	119 kDa
Uniprot ID	Q9Z2B5
Dilution Range	IHC-P=1:100-500, Flow-Cyt=2ug/Test, WB=1:500-2000, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat kidn...



Flow cytometric analysis of A431 cells I...



Flow cytometric analysis of A431 cells I...