


Product Datasheet

PHLDA3 antibody (orb50695)

Description	PHLDA3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Pleckstrin homology-like domain family A member 3 protein (43-95AA)
Target	PHLDA3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5J5
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescent analysis of PC-3 cells...