

Product Datasheet

CD1c antibody (APC) (orb507501)

Description

Mouse monoclonal antibody against CD1c conjugated to

Reactivity

Human

Conjugation

APC

Tested

FC

Applications

Immunogen

human thymocytes

Target

CD1c

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

P29017

Entrez


911

Dilution Range


Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Expiration Date


12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human CD1c
positive
CD19 p...



Flow
cytometry
surface
staining
pattern ...