



Product Datasheet

TRPM2 antibody (orb507508)


Description	Goat polyclonal antibody to TRPM2
Species/Host	Goat
Reactivity	Bovine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Target	TRPM2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Human Pancreas and Brain (Cortex) . Recommended concentration: 3.75µg/ml.. Notes: Immunofluorescence: Strong expression of the protein seen in the nuclear membranes and cytoplasm of A431 cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of Jurkat cells. Recommended
Protein Sequence	QQKQRPEQKIED
Clonality	Polyclonal
MW	171; 176; 167
Entrez	7226
Dilution Range	ELISA: 1:16000, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of paraforma...



Immunohistochemical staining of Human Co...



Immunohistochemical staining of Human Pa...