



Product Datasheet

OCIL antibody (orb507509)


Description	Goat polyclonal antibody to OCIL
Species/Host	Goat
Reactivity	Bovine, Canine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	OCIL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 30-35kDa band observed in lysates of histiocytic lymphoma cell line U937 (calculated MW of 22.3kDa according to NP_001004419.1). The observed molecular weight corresponds to earlier findings with different antibodies from other commercial sources.
Protein Sequence	LLRYKGPSDH
Clonality	Polyclonal
MW	21.8; 22.3; 17.6
Entrez	29121
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of U937 cell lysat...



Immunohistochemical staining of Human Pl...



Immunohistochemical staining of Human To...