


Product Datasheet

Human alpha-defensin 1/DEFA1 ELISA Kit (orb507514)

Description	Human alpha-defensin 1/DEFA1 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: D39-C94
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P59665
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:1000 and concentration is around 2.5 ug/ml. For plasma, dilution ratio is 1:10000 and concentration is around 30 ug/ml.
Expiration Date	Please enquire.



Human alpha-defensin/DEFA1 PicoKine ELIS...