


Product Datasheet

Mouse EFEMP1 ELISA Kit (orb507517)

Description	Mouse EFEMP1 ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q18-F493
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q8BPB5
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:400, concentration in serum and plasma is around 1600 ng/ml.
Expiration Date	Please enquire.



Mouse
EFEMP1
PicoKine
ELISA Kit
standard...