


Product Datasheet

Rat P-Cadherin-3 CDH3 ELISA Kit (orb507519)

Description	Rat P-Cadherin-3 CDH3 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E100-G647
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
Uniprot ID	F1LMI3
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 15.6 pg/ml.
Expiration Date	Please enquire.



Rat P-Cadherin PicoKine ELISA Kit standard curve