


Product Datasheet

Rat TRKB ELISA Kit (orb507521)

Description	Rat TRKB ELISA Kit
Reactivity	Rat
Range	93.6 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C32-H429
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	Q63604
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 200-700 pg/ml.
Expiration Date	Please enquire.



Rat TRKB
PicoKine
ELISA Kit
standard
cur...