


Product Datasheet

Human Dermato pontin/DPT ELISA Kit (orb507523)

Description	Human Dermato pontin/DPT ELISA Kit
Reactivity	Human
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q19-V201
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q07507
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 4-10 ng/ml.
Expiration Date	Please enquire.



Human
Dermato pontin/DPT
PicoKine ELISA K...