


Product Datasheet

Mouse PDGF-BB ELISA Kit (orb507538)

Description	Mouse PDGF-BB ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S82-T190
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	P31240
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:20 and concentration is 8216-20868 pg/ml. For plasma, dilution ratio is 1:5 and concentration is 280-9852 pg/ml.
Expiration Date	Please enquire.



Mouse
PDGF-BB
PicoKine
ELISA Kit
standar...