


Product Datasheet

Human PDGF-DD ELISA Kit (orb507541)

Description	Human PDGF-DD ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S250-R370
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q9GZP0
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:5 and concentration is around 2 ng/ml. For plasma, dilution ratio is 1:1 and concentration is 0-220 pg/ml.
Expiration Date	Please enquire.



Human
PDGF-DD
PicoKine
ELISA Kit
standar...