


Product Datasheet

Human HRG ELISA Kit (orb507545)

Description	Human HRG ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Immunogen	Expression system for standard: NS0; Immunogen sequence: V19-K525
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q6P1K1
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human HRG
PicoKine
ELISA Kit
standard
CU...