


Product Datasheet

Mouse AAT ELISA Kit (orb507547)

Description	Mouse AAT ELISA Kit
Reactivity	Mouse
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E25-K413
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q00896
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10000, concentration in serum and plasma is 30-160 ug/ml.
Expiration Date	Please enquire.



Mouse AAT
PicoKine
ELISA Kit
standard
CU...