



Product Datasheet

CD2AP Antibody (monoclonal, 5F8) (orb507557)


Description	CD2AP Antibody (monoclonal, 5F8)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E. coli-derived human CD2AP recombinant protein (Position: K253-K337).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	80 kDa
Uniprot ID	Q9Y5K6
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunocytochemistry/Immunofluorescence, 5 µg/ml Flow Cytometry, 1-3µg/1x10 ⁶ cells
Expiration Date	12 months from date of receipt.



ICC analysis of CD2AP using anti-CD2AP a...



Flow Cytometry analysis of K562 cells us...



Flow Cytometry analysis of U2OS cells us...