



Product Datasheet

SSPO antibody (orb507583)


Description	SSPO antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human SSPO. Please contact us for the exact immunogen sequence. The peptide is available as orb511020.
Target	SSPO
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	A2VEC9
Dilution Range	IHC-P: 1:50-300
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



IHC-P image of mouse stomach tissue usin...



IHC-P staining of rat cerebellum tissue ...