



Product Datasheet

DP1 antibody (orb5090)

Description	DP1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human TFDP1 (221-320/410 aa)
Target	DP1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45 kDa
Uniprot ID	Q14186
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:Kidney
(Mouse);Primary:Anti...



Western blot:Lane
1:Mouse Testis;Lane
2:...