

Product Datasheet

PNLIPRP1 antibody (orb50930)

Description	PNLIPRP1 antibody
Species/Host	Rabbit
Reactivity	Canine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Dog Inactive pancreatic lipase-related protein 1 protein (18-467AA)
Target	PNLIPRP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06857
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.