



Product Datasheet

E2F1 (phospho-Ser364) antibody (orb5095)


Description	E2F1 (phospho-Ser364) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human E2F1 around the phosphorylation site of Ser364 (MG(p-S)LR)
Target	phospho-E2F1 (Ser364)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48 kDa
Uniprot ID	Q01094
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1ug/Test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow
cytometry of
Molt4
(Black).Primary
...



ICC analysis of
HepG2 cell
using
phospho...



IHC-P analysis
of rat testis
using phospho...