



## Product Datasheet

### SLC18A2 antibody (orb51062)

|                            |                                                                            |
|----------------------------|----------------------------------------------------------------------------|
| <b>Description</b>         | SLC18A2 antibody                                                           |
| <b>Species/Host</b>        | Rabbit                                                                     |
| <b>Reactivity</b>          | Human, Mouse                                                               |
| <b>Conjugation</b>         | Unconjugated                                                               |
| <b>Tested Applications</b> | ELISA, IF, WB                                                              |
| <b>Immunogen</b>           | Recombinant Human Synaptic vesicular amine transporter protein (463-514AA) |
| <b>Target</b>              | SLC18A2                                                                    |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                |
| <b>Form/Appearance</b>     | Liquid                                                                     |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.              |
| <b>Note</b>                | For research use only                                                      |
| <b>Isotype</b>             | IgG                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                 |
| <b>Uniprot ID</b>          | <a href="#">Q05940</a>                                                     |
| <b>Dilution Range</b>      | WB:1:500-5000, IF:1:50-200                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                            |



Western blot analysis of Jurkat(lane 1),...



Immunofluorescence analysis of HepG2 cells...