


## Product Datasheet

### Chicken Ovotransferrin protein (orb510858)

|                            |                                                                         |
|----------------------------|-------------------------------------------------------------------------|
| <b>Description</b>         | Recombinant chicken Ovotransferrin protein                              |
| <b>Tested Applications</b> | ELISA, SDS-PAGE, WB                                                     |
| <b>Target</b>              | Chicken Ovotransferrin                                                  |
| <b>Preservatives</b>       | PBS, pH 7.4                                                             |
| <b>Form/Appearance</b>     | Liquid                                                                  |
| <b>Storage</b>             | Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                   |
| <b>Application notes</b>   | Tag Info: Native proteinSequence Info: Full length protein              |
| <b>Purity</b>              | > 95%(SDS-PAGE)                                                         |
| <b>Source</b>              | Purified from chicken egg white                                         |
| <b>Hazard Information</b>  | For research use only                                                   |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                          |



(Tris-Glycine gel)  
Discontinuous SDS-PAGE...