

Product Datasheet

TFEB antibody (Biotin) (orb51101)

Description	TFEB antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Transcription factor EB protein (1-476AA)
Target	TFEB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19484
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.