

Product Datasheet

SSPO peptide (orb511020)

Description	SSPO peptide
Conjugation	Unconjugated
Tested Applications	HPLC
Target	SSPO
Form/Appearance	Lyophilized
Storage	Upon reconstituted, aliquot and store at -20°C or -80°C. for long term storage, Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb507583 (CD133 antibody)
Source	Synthesis
Expiration Date	6 months from date of receipt.