



Product Datasheet

S7A13 antibody (orb511146)

Description	S7A13 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human S7A13 . Please contact us for the exact immunogen sequence.
Target	S7A13
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	52 kDa
Uniprot ID	Q8TCU3
Dilution Range	WB: 1:100-800
Expiration Date	12 months from date of receipt.



Western blot analysis of rat kidney, mouse kidney, and mouse brain tissue.



Western blot analysis of rat kidney, mouse kidney, and mouse brain tissue.