

Product Datasheet

Human Siglec-10 Protein (orb511226)

Description	The siglecs (sialic acid-binding Ig-like lectins) are a distinct subset of the Ig superfamily with a...
Reactivity	Human
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	Siglec-10
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-mFc(mlgG2a)
Note	For research use only
Application notes	This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 85.6 kDa. The protein migrates as 90-120 KDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_149121.2
Purity	95%
Source	Human Siglec-10, Mouse IgG2a Fc Tag, low endotoxin (orb511226) is expressed from human 293 cells (HEK293). It contains AA Met 17 - Thr 546 (Accession # Q96LC7-1).
MW	85.6 kDa
NCBI	NP_149121.2
Expiration Date	6 months from date of receipt.