

Product Datasheet

Cynomolgus B7-H6 / NCR3LG1 Protein (orb511255)

Description	Cynomolgus B7-H6 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Asp 25 ...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	B7-H6
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 28.4 kDa. The protein migrates as 38-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_005578557.1
Purity	95%
Source	Cynomolgus B7-H6 Protein, His Tag (orb511255) is expressed from human 293 cells (HEK293). It contains AA Asp 25 - Asp 259 (Accession # XP_005578557.1).
MW	28.4 kDa
NCBI	XP_005578557.1
Expiration Date	6 months from date of receipt.