

Product Datasheet

Human TNFSF11 / RANKL / CD254 Protein (orb511282)

Description	Human TNFSF11, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gly 64 - Asp 245 (A...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TNFSF11
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	N-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the N-terminus. The protein has a calculated MW of 46.9 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAC51762.1
Purity	90%
Source	Human TNFSF11, Fc Tag (orb511282) is expressed from human 293 cells (HEK293). It contains AA Gly 64 - Asp 245 (Accession # AAC51762.1).
MW	46.9 kDa
NCBI	AAC51762.1
Expiration Date	6 months from date of receipt.