

Product Datasheet

Human SLAMF1 / SLAM / CD150 Protein (orb511307)

Description	Self-ligand receptor of the signaling lymphocytic activation molecule (SLAM) family, also known as C...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	SLAMF1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 26.2 kDa. The protein migrates as 42-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_003028.1
Purity	95%
Source	Human SLAMF1, His Tag (orb511307) is expressed from human 293 cells (HEK293). It contains AA Ala 21 - Pro 237 (Accession # Q13291-1).
MW	26.2 kDa
NCBI	NP_003028.1
Expiration Date	6 months from date of receipt.