

Product Datasheet

Human CD52 Protein (orb511317)

Description	Human CD52, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gly 25 - Ser 36 (Acces...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD52
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 27.9 kDa. The protein migrates as 38-43 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001794.2
Purity	95%
Source	Human CD52, Fc Tag (orb511317) is expressed from human 293 cells (HEK293). It contains AA Gly 25 - Ser 36 (Accession # P31358-1).
MW	27.9 kDa
Uniprot ID	P31358
NCBI	NP_001794.2
Expiration Date	6 months from date of receipt.