

Product Datasheet

Human TSLP R Protein (orb511319)

Description	Receptor for thymic stromal lymphopoietin (TSLP). CRLF2 / TSLPR is a member of the type I cytokine r...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TSLP R
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 50.4 kDa. The protein migrates as 62-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_071431.2
Purity	95%
Source	Human TSLP R, Fc Tag (orb511319) is expressed from human 293 cells (HEK293). It contains AA Gly 25 - Lys 231 (Accession # Q9HC73-1).
MW	50.4 kDa
NCBI	NP_071431.2
Expiration Date	6 months from date of receipt.