

Product Datasheet

Mouse CD47 Protein (orb511323)

Description	Leukocyte surface antigen CD47 is also known as Antigenic surface determinant protein OA3, Integrin-...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD47
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 15.7 kDa. The protein migrates as 33-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_006521874.1
Purity	90%
Source	Mouse CD47, His Tag (orb511323) is expressed from human 293 cells (HEK293). It contains AA Gln 19 - Lys 140 (Accession # Q61735-1).
MW	15.7 kDa
NCBI	XP_006521874.1
Expiration Date	6 months from date of receipt.