

Product Datasheet

XRCC5 antibody (Biotin) (orb51149)

Description	XRCC5 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human X-ray repair cross-complementing protein 5 protein (251-455AA)
Target	XRCC5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13010
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.