

Product Datasheet

HSPA8 antibody (Biotin) (orb51157)

Description	HSPA8 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Heat shock cognate 71 kDa protein (418-597AA)
Target	HSPA8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11142
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.